

# STERILENE

**NRCS Registered - Act 5 GNR 529/237686/040/0072**

Sterilene is a pine fragranced, Glutaraldehyde, QAC blended, liquid detergent – disinfectant with a 7% total active purity. Sterilene is high foaming with anti-microbial properties for general cleaning and deodorising of ablutions, institutions, animal cages etc. Product pH @ 1 % is between 7.00 - 8.00.

## Dilutions

Use at 10 - 100 ml/L (1: 100 - 1: 10) dilution with water depending on application. Apply via spray or normal applications. Allow a 5 min contact time with a cold water rinse. Use diluted product within 7 days. Avoid contamination of foodstuffs. Do not mix product with incompatible materials, e.g. soaps and anionic detergents. Only use product on cleaned surfaces. Painted surfaces, linoleum, asphalt, rubber and plastics coverings could in some instances be adversely effected by the product and should be thoroughly rinsed with potable water within 15 min of treatment.

## Microbial Efficacy

Bactericidal – Log 5 @ 1: 100 with 5 min - *Escherichia coli*, *Staphylococcus aureus*, *Pseudomonas aeruginosa*.

Fungicidal - Log 4 – @ 1:10 with a 5 Min Contact time - *Aspergillus brasiliensis*

Fungicidal – Log 6 = @ 1: 100 with a 5 min contact time – *Candida albicans*

Sporicidal - Log 4 - @ 1:510 with a 5 min Contact time - *Bacillus subtilis spores*

## GHS Identification:

Warning

**UN 1760 – Corrosive liquid, NOS**



## PPE

Gloves, Goggles, Personal Protective Clothing, Apron, Gumboots



First Aid

SKIN	EYES	IINHALATION	ORAL
Wash off with plenty of water for 15 minutes.	Rinse cautiously with water for several minutes, remove contact lenses if safe and easy to do so, continue rinsing and get medical attention	If INHALED remove person to fresh air and keep comfortable for breathing. Do not give CPR to an unconscious person.	Do NOT induce vomiting
If skin irritation persists consult a doctor.			Neutralize with an acidic solution
	Consult a Doctor immediately	Drink Fruit Juice / Vinegar / Water	Consult a Doctor immediately

## Available in

Drums – 5Lt, 25Lt, 210Lt

Flowbin - 1000Lt